

Face Mask EarSavers



3D Printed On Demand

Instructions for Use

Made of 3D Printed PLA Plastic

Can be disinfected using:

Rubbing Alcohol, Mild Cleansers, or Soap and Water

Small EarSavers Placement Guide



Designed to be worn at the base of the head to relieve rubbing on the backs of the ears.

Medium EarSavers Placement Guide



Designed to be worn above a ponytail hairband to relieve rubbing on the backs and tops of the ears.

Designed by: Erin Shircliffe @ 3dpDesignLab